









Ropox Support Washbasin Complete

Fixed - Basic - Manual USER MANUAL MOUNTING INSTRUCTIONS

Keep this manual with the product at all times

Contents:

1.	INTRODUCTION	
2.	CONFORMITY WITH EU DIRECTIVES	3
3.	APPLICATION	3
	3.1 Operation	4
4.	TECHNICAL DATA	5
5.	DIMENSION SKETCHES	6
6.	INSTALLATION INTERFACES	7
	6.1 Connection adapters	11
7.	MOUNTING INSTRUCTIONS	12
	 7.1 Support Washbasin, Fixed (40-44000/01) 7.2 Support Washbasin with Basic frame (40-4002/03) 7.3 Mounting of covers for 40-44002/03 7.4 Support Washbasin with Manual frame (40-44004/05) 7.5 Mounting of covers for 40-44004/05 	
8.	LIST OF COMPONENTS	33
9.	PERFORMANCE TEST	37
10.	SAFETY IN USE	37
11.	SERVICE AND MAINTENANCE	38
	 11.1 Cleaning	38 38
12.	FAULT FINDING	39
13.	LIST OF SPARE PARTS	40
14.	OPTIONAL EXTRAS	40
15.	LABEL	41
16	COMPLAINT	42

1. Introduction

Owing to its unique design, Ropox Support Washbasin provides additional support and easier access to the washbasin for wheelchair users and walking-impaired persons. The manual height adjustment facilitates frequent adjustment of the washbasin to suit several users. The height adjustment range of the adjustable models is 30cm.

The ample space for knees and legs makes the washbasin ideal for wheelchair users, who obtain optimum mobility and a perfect working position.

The integrated grab handles facilitate access and allows the user to get close to the washbasin.



This manual must ALWAYS accompany the product and be read by and available to the user

Correct use, operation and inspection are decisive factors for effective and safe performance.

2. Conformity with EU Directives

This product has CE marking and is thus in accordance with the basic requirements concerning functionality and safety according to the current Machinery Directive.

All models of Ropox Support Washbasin are manufactured in accordance with the current requirements of Directive No. 93/42/EEC concerning Medical Devices and belong to the lowest risk class 1

3. Application


Ropox Support Washbasin Basic and Manual may be adjusted in height for optimum working height for user and carer. Ropox Support Washbasin must NOT be used as a lifting table or person lifter.

The product is designed for indoor use under normal room temperature and humidity as specified in Section 4.

3.1 Operation


Basic height adjustment (applying to 40-44002/03)

The height of the washbasin may be adjusted by loosening the two handles. Note that the washbasin must be supported while loosening the handles in order to prevent it from falling to the bottom position. Covers, if any, must be removed before adjustment of these models.



Manual height adjustment (applying to 40-44004/05):

The height of the washbasin may be easily adjusted by means of the detachable hand crank.



Technical data

Product name:

tom number: Hom number Name Illustration

Ropox Support Washbasin

Item number: Name Illustration

40-44000 Fixed (left)

40-44001 Fixed (right)

40-44002 Basic height adjustment (left)

40-44003 Basic height adjustment (right)

40-44004 Manual height adjustment (left)

40-44005 Manual height adjustment (right)

Adjustment range: 30 cm (71-101cm)

Mixer tap: Oras Saga, mixer tap with temperature pre-setting facility.

Alternative mixer taps available.
Contact Ropox for more information.

Materials: Welded steel pipes St. 37

Aluminium alloy 6060-T6

Various plastic components

Washbasin: Synthetic/marble composite Cover: Vacuum-shaped ABS

Surface treatment: Powder coating, RAL 9010 matte white, chromite, gel coat

Max. load: Tested in accordance with DS/ISO 17966:2016, max user weight is

400kg/881lb

Ambient temperature

range:

-5 °C to +40 °C Transport and storage

+5 °C to +40 °C Operation

Relative humidity: 20% to 80% - non condensing

Complaint: See Complaint page 42

Manufacturer: Ropox A/S, DK-4700 Næstved, Phone: +45 55 75 05 00


e-mail: info@ropox.dk, www.ropox.dk

Intended use


The support washbasin is especially suitable for walking-impaired because you can use the integrated, ergonomic grab rails as support and as help when getting up or sitting down. The washbasin is available with fixed height and with height adjustment. It is possible to choose between manual or electrical height adjustment. This makes it able to suit the individual user who will be using the washbasin.

5. Dimension sketches


Fixed 40-44000/01




Basic height adjustment 40-44002/03



Manual height adjustment 40-44004/05



Accessories – covers 40-44021/22



Page 7

© ROPOX 2017

6. Installation interfaces

Fixed (40-44000 and 40-44001)

CONSTRUCTION

The wall must be made of a material suited for screw mounting and approved for damp rooms. The hatched area is the minimum area for reinforcement of the wall. Reinforcement from floor to ceiling is advisable.

Load-tested in Accordance to DS/ISO 17966, the maximum user weight has been fixed at 400kg (881 lbs).

Max. tensile load of upper screws 150kg/screw (330 lbs)

The washbasin is supplied with a set of 6 screws and rawlplugs suitable for mounting on a concrete wall. The fitter should always consider the material, condition and strength of the wall and use screws and rawlplugs suitable for the specific wall type.

PLUMBING

The drain is \emptyset 32 (1 ½") and ends in an \emptyset 32x88.5° drain elbow.

The maximum protrusion of drain connection and elbow from the wall is 90mm (3 $\frac{1}{2}$ ") to allow for the covers. The water supply ends in 2 (two) 90° shut-off valves with $\frac{1}{2}$ " outside thread pointing downward (not included). The maximum protrusion of the shut-off valves from the wall is 90mm (3 $\frac{1}{2}$ ") to allow for the covers.


If a location of the water supply and drain different from the recommended positioning is required or if there is a floor drain, the Support Washbasin will have no covers.


For purchase of adapters for connection to various standard installations, see "Connection Adapters – Other Washbasins".


NOTE:

All dimensions are based on a washbasin height of 710mm (28") above floor level. If a change of this height is required, the recommended dimensions must be changed accordingly.

The illustrations show the left-hand model 40-44000







Basic

(40-44002 and 40-44003)

CONSTRUCTION

The wall must be made of a material suited for screw mounting and approved for damp rooms.

The hatched area is the minimum area for reinforcement of the wall.

Reinforcement from floor to ceiling is advisable.

Load-tested in Accordance to DS/ISO 17966, the maximum user weight has been fixed at 400kg (881 lbs).

Max. tensile load of upper screws 85kg/screw (187 lbs)

The washbasin is supplied with a set of 6 screws and rawlplugs suitable for mounting on a concrete wall. The fitter should always consider the material, condition and strength of the wall and use screws and rawlplugs suitable for the specific wall type.

PLUMBING

The drain is \emptyset 32 (1 $\frac{1}{4}$ ") and ends in an \emptyset 32x88.5° drain elbow.


The maximum protrusion of drain connection and elbow from the wall is 90mm (3 $\frac{1}{2}$ ") to allow for the covers. The water supply ends in 2 (two) 90° shut-off valves with $\frac{1}{2}$ " outside thread pointing downward (not included). The maximum protrusion of the shut-off valves from the wall is 90mm (3 $\frac{1}{2}$ ") to allow for the covers.


It is important that all hoses for water connection and drain are flexible to allow the frame to move freely within the height adjustment range of 710-1010mm (28" – 39 8/10") (including washbasin).


If a location of water supply and drain different from the recommended positioning is required or if there is a floor drain, the Support Washbasin will have no covers. For purchase of adapters for connection to various standard installations, see "Connection Adapters – Other Washbasins".

NOTE:

The illustrations show the left-hand model 40-44002







Manual (40-44004 and 40-44005)

CONSTRUCTION

The wall must be made of a material suited for screw mounting and approved for damp rooms.

The hatched area is the minimum area for reinforcement of the wall.

Reinforcement from floor to ceiling is advisable.

Load-tested in Accordance to DS/ISO 17966, the maximum user weight has been fixed at 400kg (881 lbs).

Max tensile load of upper screws 71kg/screw (157 lbs)

The washbasin is supplied with a set of 8 screws and rawlplugs suitable for mounting on a concrete wall. The fitter should always consider the material, condition and strength of the wall and use screws and rawlplugs suitable for the specific wall type.

PLUMBING

The drain is \emptyset 32 (1 $\frac{1}{4}$ ") and ends in an \emptyset 32x88.5° drain elbow.


The maximum protrusion of drain connection of elbow from the wall is 90mm (3 ½") to allow for the covers. The water supply ends in 2 (two) 90° shut-off valves pointing downward (not included). The maximum protrusion of the shut-off valves from the wall is 90mm (3 ½") to allow for the covers


It is important that all hoses for water connection and drain are flexible to allow the frame to move freely within the height adjustment range of 710-1010mm (28" – 39 8/10") (including washbasin). If a location of water supply and drain different from the recommended positioning is required and if there is a floor drain, the


For purchase of adapters for connection to various standard installations, see "Connection Adapters – Other Washbasins".

Support Washbasin will have no covers.

NOTE: The illustrations show the left-hand model 40-44004







6.1 Connection adapters
The highlighted components can be bought separately depending on the type of building installation requirements.

Water supply

Hoses from basin	Adapters	Installation in room (not included)
	No adapter necessary	1/2 " Thread
1/2"		
	1/2" -> 3/8" Item no. 97001670	3/8 " Thread
	1/2" -> 15mm	15mm
	Item no. 97001666	Pipe

Water hose extensions

Water 1103c exterisions		
Hose from basin	Coupling piece	Extra hose
		1/2" x 1/2", 0,5m, steel
		Item no. 97001123
	1/2" -> 1/2" Item no. 97001662	1/2" x 1/2", 1m, plastic Item no. 97001120

Drain


Diaiii		
Drain pipe from basin	Adapters	Installation in room (not included)
Ø32	No adapter necessary	Ø32
	Ø32 -> Ø40 Item no. 97001660	Ø40
	Ø32 -> Ø50 Item no. 97001661	Ø50

Drain pipe extensions


Drain pipe from basin Coupling piece Extra drain pipe		
Drain pipe nom basin	Coupling piece	Extra drain pipe
	Ø32 Flex hose, 27-62cm Item no. 97001162	The state of the s
Ø32		
	Ø32-Ø32 Item no. 97001662	Ø32 Flex hose, 40-110cm Item no. 97001161


7. Mounting instructions


7.1 Support Washbasin, Fixed (40-44000/01)




Mounting of covers: See page 17-20










Fine adjustment of wall fitting – use legs on washbasin as a guide














7.3 Mounting of covers for 40-44002/03





















7.4 Support Washbasin with Manual frame (40-44004/05)





















Tip: Move washbasin to top position to facilitate further mounting







7.5 Mounting of covers for 40-44004/05







Page 29














8. List of components

Ropox Support Washbasin, Fixed, left, 40-44000:

40*44010-064: Top plate with washbasin, left-oriented 1 (one)

Provided with:

40*44010-063: Fixed wall bracket, long 1 (one)

40*44010-062: Fixed wall bracket, short 1 (one)

40*44010-048: Frame Basic side bracket 1 (one)

40-44021: Cover 1 (one) set

Ropox Support Washbasin, Fixed, right, 40-44001:

40*44010-070: Top plate with washbasin, right-oriented 1 (one)

Provided with:

40-44021:

40*44010-063: Fixed wall bracket, long 1 (one)

40*44010-062: Fixed wall bracket, short 1 (one)

40*44010-048: Frame Basic side bracket 1 (one)

40-44022: Cover 1 (one) set

Ropox Support Washbasin, Basic height adjustment, left, 40-44002:

40*44010-064: Top plate with washbasin, left-oriented 1 (one)

Basic height-adjustment brackets 1 (one) set

 Incl.
 40*44010-025
 Long bracket

 40*44010-026
 Short bracket

 40*44010-020
 Inside legs

 97800-022
 Threaded plate

 97001631
 Adjustable handles

 40*44010-048
 Frame Basic side bracket

Cover

Ropox Support Washbasin, Basic height adjustment, right, 40-44003:

40*44010-070: Top plate with washbasin, right-oriented 1 (one)



Basic height-adjustment brackets

Incl. 40*44010-025 40*44010-026 40*44010-020 97800-022 97001631

40*44010-048

Long bracket Short bracket Inside legs Threaded plate Adjustable handles

Frame Basic side bracket

40-44022: Cover



1 (one) set



Ropox Support Washbasin, Manual height adjustment, left, 40-44004:

40*44010-064: Top plate with washbasin, left oriented 1 (one)



Provided with:

Top frame, left oriented

Incl.

40*44010-042 40*44010-057

Frame Manual support, left Hand crank assembly left

1 (one)



30-67600:

Flexi legs, Manual

Incl.

Distance blocks

Cover

2 (two) sets



40*44010-038:

Hexagonal shaft

1 (one)



20*60320-548:

Hand crank

1 (one)



40-44027:

Cover

1 (one) set



Ropox Support Washbasin, Manual height adjustment, right, 40-44015:

40*44010-070: Top plate with washbasin, right-oriented 1 (one)



Provided with:

Top frame, right-oriented

1 (one)

Incl. 40*44010-077 40*44010-082

Frame Manual support, right Hand crank assembly, right



30-67600: Flexi legs, Manual

Cover

Distance blocks



40*44010-038:

Hexagonal shaft

1 (one)

20*60320-548:

Hand crank



40-44028:

Cover

1 (one) set



8.1 Common parts

40-44020	Water and drain kit	
Incl. 97001601	Siphon	1 (one)
97001162	Flexible drain hose ø32/32 pipe stub	1 (one)
97001035	ø32x88.5° angle	1 (one)
97001120	Connecting hose, ½" x ½", 100 cm:	2 (two)
40*44010-098	Drain relief	1 (one)
40-44040	Mixer tap, Oras Saga	1 (one)

9. Performance test

After completed installation and before use, all functions of the Ropox Support Washbasin must be tested. After that, a performance test must be carried out at least once a year by competent personnel.

- 1. Check that the mounting instructions have been observed.
- 2. Check that all bolts and screws have been tightened.
- 3. Check that all drains and hoses have been correctly connected.
- 4. Ensure that there is no load on the frame.
- 5. Ensure that there are no objects preventing the frame from moving freely within the height adjustment range.

When these tests have been carried out successfully, Ropox Support Washbasin is ready for use. See Section 10 "Safety in use".

10. Safety in use

- Ropox Support Washbasin must only be adjusted in height by persons who have read and understood these instructions.
- Some models of Ropox Support Washbasin are height-adjustable and must not be used as a lifting table or person lifter.
- Always use the washbasin in a manner excluding the risk of damage to persons or property. The person operating the washbasin is responsible for avoiding damage or injury.
- If the washbasin is used in publicly accessible locations where children or persons with reduced observation capacity may get near it, the person operating the washbasin must pay attention to those present in order to prevent dangerous situations.
- Make sure that there is free space above and below the washbasin to allow height adjustment.
- Do not overload Ropox Support Washbasin and make sure that the load distribution is correct.
- > Do not operate the washbasin in case of defects or damage.
- Modifications, which may influence the performance or design of Ropox Support Washbasin, are not allowed.
- Installation, service and repairs must be carried out by competent personnel.
- During inspection, service or repairs any load must be removed from Ropox Support Washbasin.

Service and maintenance

11.1 Cleaning

During use, the product may get dirty. It is very important to clean the product as described in these instructions.

Wipe Ropox Support Washbasin with a wrung cloth with lukewarm water and a mild detergent. Do not use enchants, abrasives or grinding cloths, brushes and sponges.

11.2 Cleaning of tabletop and washbasin

Ropox Support Washbasin is made of a composite material of marble and plastic granulate with a gel coat surface. This gives a hard and very hardwearing finish, which is easy to clean. Wash the surfaces daily with water and dry with a cloth. Do not use detergents containing hard abrasives – e.g. steel wool, scouring powder or sulphuric acid.

If the water has high calcium content, calcareous deposits will often occur. These deposits should be removed at least once a week by using a detergent containing citric acid, acetic acid or phosphoric acid. Leave the detergent to work for 5 to 10 minutes depending on the thickness of the deposits and clean thoroughly with a soft scouring sponge.

Disinfection of the basin may be done with household spirits.

11.3 Inspection



Inspections, service and repairs must be carried out by competent personnel

The frame of Ropox Support Washbasin is maintenance-free. Moving parts have been lubricated for life. For reasons of safety and reliability, we recommend inspection of the frame once a year.

Annual inspection:

- Check that all bolts and screw have been securely tightened.
- Check that the frame moves freely and without problems from bottom to top position.
- Check that water hoses and drain hoses are tight and functional.

After each inspection, complete the service schedule. See Section 11.5.



Only use Ropox original spare parts for replacement. If other spare parts are used, the whole guarantee may become void.

11.4 Service schedule

Service and maintenance Serial No.	Service and maintenance Serial No.
Date:	Date:
Signature:	Signature:
Comments:	Comments:
Service and maintenance Serial No.	Service and maintenance Serial No.
Date:	Date:
Signature:	Signature:
Comments:	Comments:

12. Fault finding

a) The frame seems loose or unstable

- 1. The assembly screws of the frame have not been securely tightened.
- 2. Retighten all screws of the frame according to the mounting instructions.

b) The washbasin cannot be adjusted in height / adjustment is very difficult

1. Check that there are no jammed objects preventing adjustment.

c) It takes too long for the water to leave the washbasin

- 1. Check the siphon for clogging and clean it, if required.
- 2. Check that the flexible hose on the wall has not been positioned so high that it acts as an extra siphon. Cf. recommendation, Section 6.

13. List of spare parts


Handle M6 Item No. 96000177 Applies to 40-44012/13 only	
Hand crank Item No. 20-60320-548 Applies to 40-44014/15 only	

14. Optional extras

Ropox overflow plug Item No. 40-40935	

15. Label

Each unit has a label of its own as illustrated. This label contains important information to be stated in connection with inquiries to the supplier. For positioning of the label, see drawing of spare parts.



	This product is CE-marked in accordance with:
CE	European Medical Device Directive 93/42/EEC, including amendments incorporated in Directive 2007/47/EEC.
	European Medical Device Regulation (EU)2017/745
	Council Directive 2011/65/EU, RoHS
	Manufacturer name and address
M	Date of production
REF	Stock number
SN	Serial number
Ţ <u>i</u>	Consult manual before use
<u> </u>	Consult manual for important safety related information, warnings and safety precautions.
Z	Do not dispose as unsorted municipal waste. Product must be returned to a designated recycling station.
5°C 40°C	Operating temperatures may not exceed 5-40C°

16. Complaint

See General Terms of Sale and Delivery on www.ropox.dk



Ringstedgade 221 DK – 4700 Næstved Tel.: +45 55 75 05 00 Fax.: +45 55 75 05 50 E-mail: <u>info@ropox.dk</u> <u>www.ropox.dk</u>