

ROPOX VARIO CHANGING BED

SPACESAVING CHANGING BED

The Ropox Vario changing bed is mounted on the wall.

The height is adjustable from 30-100cm using the hand control with spiral line.

When the changing bed is not in use, the tabletop can be folded up and locked in a vertical position so that it takes up very little space in the room.

Folding up and down the table top takes place in a muted movement.

Ropox Vario changing bed is available in five models incl. and without a folding side bed guard and mattress adapted to the changing table.

The mattress has a cleaning-friendly surface that makes it easy to clean after use.







ROPOX VARIO CHANGING BED

CHANGING BED

ITEM NO.	DESCRIPTION
40-30604	Vario without bed guard rail 120x70cm
40-30606	Vario without bed guard 140x70cm
40-30608	Vario without bed guard 160x70cm
40-30610	Vario without bed guard 180x70cm
40-30611	Vario without bed guard 190x70cm


ITEM NO.	DESCRIPTION
40-30804	Vario with bed guard 120x70cm
40-30806	Vario with bed guard 140x70cm
40-30808	Vario with bed guard 160x70cm
40-30810	Vario with bed guard 180x70cm
40-30811	Vario with bed guard 190x70cm

ACCESSORY

ITEM NO.	DESCRIPTION
40-25782	Cushion to the bed guard length 120/140cm
40-25783	Cushion to the bed guard length 160/180/190cm

ITEM NO.	DESCRIPTION
40-25760	Matress for 120x70cm
40-25795	Matress for 140x70cm
40-25796	Matress for 160x70cm
40-25764	Matress for 180x70cm
40-25765	Matress for 190x70cm

MÅLSKITSER





The items are produced according to the following Standard: DS/EN 12182:2012.