# Designed to make your life easier"

# ROPOX ERGO TABLE

## A FUNCTIONAL TABLE THAT SATISFIES YOUR EVERY NEED

Ergo Tables have a broad range of adjustment options, which allows you to adjust the table according to your needs. The variety also makes them highly suitable as tables for training, teaching, work, and play.

You can adjust the table height manually, and since the table has no cross bars, wheelchair users can easily get close access to the table.

The tiltable table top makes it easy to reach the worktop, and it also shortens the visual distance to the documents on table. This is a great benefit for weak-sighted, users with reduced arm function and people who have difficulties leaning over a table.

If required, you can use arm supports to relieve your forearms and provide a slight body support for a more balanced sitting position.




#### Max. load capacity:

# ROPOX ERGO TABLE

Height: Depth: Width: Inner Width: Frame Colour: 100 kg (on ruler: about 8 kg) 56-90 cm 60 cm 90 or 120 cm 73.6 or 104.6 cm Silver Grey

### TECHNICAL INFORMATION



All tables are manually height adjustable from 56-90 cm.

An additional height extension can be used for an adjustment interval from 66-100 cm.

The tabletop can be tilted with stepped angles up to 54° (B to D variant).

### ACCESSORIES

Arm support incl. bag and hook Item: 20-60520



Additional height extension, 10 cm Item: 20-70503



Castors with brake without height increase Item: 20-70502



Arm support lock-able incl. hook and bag Item: 20-60517



Document/PC holder incl. mounting plate Item: 20-70504



4 Brake wheels, increase height by 9 cm Item: 20-50300



Supporting ruler Item: 20-70509/90 cm (B) Item: 20-70510/120 cm (B) Item: 20-70513/90 cm (C+D) Item: 20-70514/120 cm (C+D)



Pen holder Item: 20-70506



4 Brake wheels, increase height by 8 cm Item: 20-50340



