

”When you want  
to spare your efforts”



## ROPOX SHOWER SEAT

### ALLOWS YOU TO TAKE A REST IN THE SHOWER

A shower seat enables you to sit safely and comfortably while showering.

The seat is foldable, therefore it takes up very little space when it is not in use. It has a very soft padding and is made in a hard-wearing, insulating material which does not feel cold on the skin.

The seat has a large surface area and a solid grabbing edge so that you can have a comfortable rest in safe surroundings.

You can easily adjust the seat's height to meet individual needs, and if you would like more support, you can choose a backrest to lean on either with or without arm supports.







# ROPOX SHOWER SEAT

Height  
adjustment: 12,5 cm

Max. user weight: Tested acc. to DS/ISO 17966:2016  
(see max. user weight below)

## VARIANT/MODEL


SHOWER SEATS	
	
With leg 200 kg Item no.: 40-43011	Without leg 150 kg Item no.: 40-43016

BACKRESTS	
	
With arm support Item no: 40-43025	Without arm support Item no.: 40-43030

The shower seat and the arm supports are foldable which provides more space in the bathroom.

## ACCESSORIES

Height adjustable unit for shower  
seat and backrest  
Item no.: 40-43043



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS/ISO 17966:2016.

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

**ROPOX**<sup>TM</sup>  
BETTER WAYS TO BETTER DAYS  
— A PART OF AddLife

No.1962308-0221