

”A bathroom aid that extends independence ”

## ROPOX TOILET LIFTER

### MAKING IT EASIER TO USE THE TOILET INDEPENDENTLY

The ROPOX Toilet Lifter is a bathroom aid that gently helps a user on and off the toilet. It is a small unit with high capacity and a discreet design that is easily customised to individual user needs.

It comes in a manual and an electrical version with height adjustment ranges of 22 and 25 cm respectively. The electric model is adjusted with a wireless hand control. The manual model is operated with a handle which can be demounted if desired.

Height-adjustable Toilet Support Arms can be mounted on the side of the unit. They provide plenty of space and good access for transfers. They can be installed or removed one at a time if user needs require it.

The smooth surface with rounded edges makes it easy to clean and helps to ensure a high level of hygiene.



**ROPOX™**

BETTER WAYS TO BETTER DAYS



- A PART OF AddLife®

# ROPOX TOILET LIFTER

Height adjustment: Electric 25 cm (from 40-65 cm)  
Manual 22 cm (from 40-62 cm)

Maximum load: Max. user load cf.  
DS / ISO 17966:2016

## VARIANT/MODEL

TOILET LIFTER	
	
Electric Item no.: 40-45027	Manual Item no.: 40-45026


Toilet lifter unit in bottom position: Model Electric/Manual 115,4x 52,3 x 20,5 cm.

Geberit "top of the line" cisterne components.

Dual flushing system that can be adjusted for "green building" requirements with 2 and 4,5 liter flush.

## ACCESSORIES


Height adjustable  
Toilet support arms, Wave  
Item no.: 40-45013 (66 cm)  
Item no.: 40-45014 (88 cm)



Toilet roll holder  
Item no.: 40-40910



Height adjustable  
Toilet support arms, Straight  
Item no.: 40-45011 (66 cm)  
Item no.: 40-45012 (88 cm)




Standing/Sitting attachment  
Item no.: 40-40925



Wall-hung toilet, Short model  
Extra dirt-repellent surface  
V-nr: 40-44070  
(depth 56 cm)



Long backrest  
Item no.: 40-45100



Wall-hung toilet, Long model  
Extra dirt-repellent surface  
V-nr: 40-44071  
(depth 70 cm)



Short backrest  
Item no.: 40-45101



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS / ISO 17966:2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

**ROPOX**<sup>TM</sup>  
BETTER WAYS TO BETTER DAYS  
- A PART OF AddLife

No.1962309-0221