

”A bathtub everyone  
can use”



## ROPOX BATHTUB

### EASY AND GENTLE BATHING

When you need extra support from a carer to take a bath, a bathtub that can be adjusted in height makes it easy and comfortable to get in and out of the bath.

It is easy to adjust the tub to the optimum working height for the helpers. The helpers can also use a mobile lift for helping a person to the bathtub as there is the necessary room under the tub for the lifting device. This provides safety and satisfaction for both carer and user.

The height is adjusted with the hand control and getting in and out of the bath matches the user's physical capabilities, this also provides the optimum working position for the helper.

A shower bed cover can be fastened to the top of the tub adding to the comfort when bathing.

The bathtub is made in a simple design making it easy to clean and giving a high level of hygiene.

The Ropox frame is extremely strong and has been tested with much more weight, than the max user weight indicated, thereby ensuring the absolute highest level of safety for the user.



**ROPOX™**




BETTER WAYS TO BETTER DAYS

- A PART OF AddLife®

# ROPOX BATHTUB


Height adjustment: 30 cm (61-91 cm) above floor  
 Maximum load: 250 kg (incl. user and water)  
 cf. DS/ISO 17966:2016 350

## VARIANT/MODEL

BATHTUB - 160X70 CM 	BATHTUB - 170X75 CM 	BATHTUB - 180X80 CM 
		
Inside measurements 142x57 cm Item no.: 40-14156	Inside measurements 152x62 cm Item no.: 40-14157	Inside measurements 164x67 cm Item no.: 40-14158

The tub is 40 cm deep on all models.

Bathtub incl. hand shower, mixer tap and flexible plumbing kit.



The bathtub mounted with a foldable shower bed (model chosen based on location of mixer tap) and bed guard.

The bathtub is easy to heightadjust with the electrical hand control which is included.

## ACCESSORIES

Shower bed for mounting on the bathtub

Item no.: 40-14226 (right)/40-14236 (left) for 160 cm

Item no.: 40-14227 (right)/40-14237 (left) for 170 cm

Item no.: 40-14228 (right)/40-14238 (left) for 180 cm

Support for hand shower

Item no.: 40-14120

To be ordered with the tub

Bed guard for bathtub

Item no.: 40-14186 for 160 cm

Item no.: 40-14187 for 170 cm


Item no.: 40-14188 for 180 cm

Brake wheels incl. fittings (4 pcs.)

Item no.: 40-14111 for 160 cm

Item no.: 40-14112 for 170 cm

Item no.: 40-14113 for 180 cm




Safety stop (one side)

Item no.: 40-14206 for 160 cm

Item no.: 40-14207 for 170 cm

Item no.: 40-14208 for 180 cm





Safety stop (two sides)

Item no.: 40-14216 for 160 cm

Item no.: 40-14217 for 170 cm

Item no.: 40-14218 for 180 cm



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS / ISO 17966:2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark – Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

**ROPOX™**

BETTER WAYS TO BETTER DAYS

– A PART OF AddLife®