

”Become more self-sufficient  
with a little extra support”



## ROPOX GRAB AND SHOWER RAILS

### SAFELY MANOEUVRE WITH GRAB RAILS

Grab and shower rails are robust and stable aids for the independent user in situations where they need a little extra support.


The grab rails are available in a number of lengths and shapes enabling you to design a bathroom with the right aids at the right place for the individual user.

Due to the design and size of the grab rails, you will have a firm and stable grip. They are easy to clean thus keeping a high hygiene level.


The grab rails for adjusting the height of the hand shower can be used with just one hand making it very easy to use.




## VARIANT/MODEL

GRAB RAILS				
				
45 cm	70 cm	80 cm	100 cm	80 x 60 cm
Item no.: 40-44608 white Item no.: 40-44618 grey	Item no.: 40-44609 white Item no.: 40-44619 grey	Item no.: 40-44610 white Item no.: 40-44620 grey	Item no.: 40-44614 white Item no.: 40-44624 grey	Item no.: 40-44615 white Item no.: 40-44625 grey

SHOWER RAILS WITH SHOWER HOLDER	
	
100 cm	120 x 45 cm
Item no.: 40-44611 white Item no.: 40-44621 grey	Item no.: 40-44612 white Item no.: 40-44622 grey


  


55 x 10 x 100 cm
Item no.: 40-44613 white Item no.: 40-44623 grey

The modular grab and shower rails have a special oval shape which enables a natural firm grip. They are produced in aluminium and polymer plastic, with a non-slip, anti-bacterial surface.

## ACCESSORY

Shower shelf  
Item no.: 40-44130



Adjustment pin for shower holder Increases the distance between adjustment knob and shower holder by 35cm  
Item no.: 40-44600



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.