### "Making it safe and comfortable to take a shower"



## ROPOX SHOWER/CHANGING BED

### SAFE AND EASY TO USE

The shower/changing bed gives the user optimum comfort and safety, and the helper gets the best possible working position. This provides security and satisfaction for both parties.

The bed is made of double-sided slats. One side of the slat curves slightly downwards ensuring a safe position for the person bathing. The other side curves slightly upwards making it easy for a helper to turn the user on the side when using the bed for changing clothes.

It is easy to adjust the height of the bed with the hand control, and the bed can quickly be adjusted to the best position, e.g. for a wheelchair so transferring a user onto the bed is smooth and easy.

The shower/changing bed is cleaning-friendly, and can be folded up when not in use to save space in the bathroom.







# **ROPOX SHOWER/CHANGING BED**

Height adjustment: 70 cm (30-100 cm) above floor

Max. user weight: Testet acc. to DS/ISO 17966:2016 Max. user weight 220 kg

### VARIANT/MODEL



- The bed can be folded up to save space in the bathroom.
- The shower/changing bed includes an electrical hand control for easy height adjustment.
- The lying area is 18 cm shorter than the length of the shower/changing bed. •

#### **ACCESSORIES**

#### Bed guard:

Item no.: 40-25034 for shower/changing bed 146 cm Item no.: 40-25038 Item no.: 40-25036 for shower/changing bed 178 cm Item no.: 40-25035 for shower/changing bed 210 cm

End guard:

Back support: Item no.: 40-25031

Headrest: Item no.: 40-25090









NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

