

”With the right aids,  
daily life becomes easier”



## ROPOX SUPPORTLINE

### A ROBUST WASHBASIN THAT CAN BE ADJUSTED TO MATCH YOUR NEEDS

A support washbasin is especially suitable for walking-impaired because you can use the integrated, ergonomic grab rails as support and as help when getting up or sitting down.

You can adjust the height of the washbasin to suit the physical capabilities of the individual user, and it is easy to slide the unit sideways allowing you to place the basin where you need it the most.

It is easy to keep the washbasin clean due to the cleaning-friendly design, and even though the basin is very robust and stable, it still has a simple and elegant look.

The solid, integrated grab rails ensure safety to lean on and use as a support when you move around in the room.













**ROPOX™**  
BETTER WAYS TO BETTER DAYS  
— A PART OF AddLife®

# ROPOX SUPPORTLINE

Height adjustment: Model Basic / Manual 71-101 cm  
Model Electric 71-96 cm

Max. user weight: Tested acc. DS/ISO 17966:2016  
(see max. user weight below)


## VARIANT/MODEL

COMPLETE ELECTRIC 		COMPLETE MANUAL 		COMPLETE BASIC 	
					
Left / 520 kg Item no.: 40-44006	Right / 520 kg Item no.: 40-44007	Left / 400 kg Item no.: 40-44004	Right / 400 kg Item no.: 40-44005	Left / 400 kg Item no.: 40-44002	Right / 400 kg Item no.: 40-44003
ELEKTRIC HEIGHT ADJUSTMENT 		MANUAL HEIGHT ADJUSTMENT 		BASIC - PERIODIC ADJUSTMENT 	
					
Left / 520 kg Item no.: 40-44016	Right / 520 kg Item no.: 40-44017	Left / 400 kg Item no.: 40-44014	Right / 400 kg Item no.: 40-44015	Left / 400 kg Item no.: 40-44012	Right / 400 kg Item no.: 40-44013
COMPLETE FIXED 		FIXED - FIXED HØJDE 			
					
Left / 400 kg Item no.: 40-44000	Right / 400 kg Item no.: 40-44001	Left / 400 kg Item no.: 40-44010	Right / 400 kg Item no.: 40-44011		


## ACCESSORIES

### Cover

Item no.: 40-44023 (left)+40-44024 (right) (Electric)  
Item no.: 40-44027 (left)+40-44028 (right) (Manual)  
Item no.: 40-44021 (left)+40-44022 (right) (Basic & Fixed)




System for sideways sliding (35 cm)  
Item no.: 40-44030 (Manual & Electric)




### Plumbing kit (next to unit)

Item no.: 40-44034 (Electric & Manual)



### Plumbing kit (within unit)

Item no.: 40-44031 (Manual)



### Plumbing kit

Item no.: 40-44020+31 (Manual)  
Item no.: 40-44033 (Basic & Fixed)



### Fixed press pad

Item no.: 40-44026



### Plug with overflow

Item no.: 40-40935



### One-handle mixer tap

Item no.: 40-44040



### Extended lever for one-handle mixer tap

Item no.: 40-44041



### Touchless mixer tap

Item no.: 40-44045



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS / ISO 17966:2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

# ROPOX™

BETTER WAYS TO BETTER DAYS

- A PART OF AddLife®