

ROPOX SWINGLINE

A FLEXIBLE WASHBASIN THAT MAKES IT EASY TO WASH HANDS

The swing washbasin is recommended for both seated and standing users. You can adjust the height and rotate it 180°, and you can pull it to the user.

Walking-impaired users are particularly content with this washbasin because it can be pulled to the toilet. This means that they do not have to get up to use it, but instead, they can safely wash their hands while seated on the toilet.

The swing washbasin is also very practical for wheel-chair users. You can move the basin when required, e.g. push the washbasin to the corner when you need extra room for turning the wheelchair.

The swing washbasins with a lifting unit can easily be adjusted in height when needed. The simple design without any dirt traps makes it easy to clean.








ROPOX SWINGLINE


Height adjustment: 15 cm without Dock-In 20 cm with Dock-In


Max. user weight: Tested acc.to DS/ISO 17966:2016

max. user weight 245 kg in and 80 kg without Dock-In

VARIANT/MODEL









ACCESSORIES

Accessory kit without wall bracket/with wall bracket Item no.: 40-41145/40-41150



Plug with overflow Trap Accessory tray Item no.: 40-41146 Item no.: 40-40935





Extended lever for Item no.: 40-41170 one-handle mixer tap Item no.: 40-44041



Touchless mixer tap Item no.: 40-44045



Please note that all parts are suitable for dish washer/Autoclave.

NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.