

”More independence
with your personal care”



ROPOX TOILET SUPPORT ARMS

GET SOMETHING TO LEAN ON

Toilet support arms are a great aid if you need to get on and off the toilet safely. Their ergonomic shapes and rubber surface provide a comfortable grip, and the cleaning-friendly surfaces optimize hygiene.

Because the Wave models are designed with two heights and an anti-slip surface, you have a firm grip from both standing and seated position.

The support arms can be adjusted in height to suit the individual user's needs.

If you need extra room in the bathroom, the support arms can be folded up to vertical position.



ROPOX™

BETTER WAYS TO BETTER DAYS








- A PART OF AddLife®

ROPOX TOILET SUPPORT ARMS

Height adjustment: 20 cm (with height adjustabel unit)

Max. user weight: Tested acc. to DS/ISO 17966:2016 (see max. user weight below)

VARIANT/MODEL

WAVE - 90 CM		WAVE - 76 CM	
			
250 kg Item no.: 40-40421 grey Item no.: 40-40426 red	330 kg Item no.: 40-40431 grey Item no.: 40-40436 red	250 kg Item no.: 40-40420 grey Item no.: 40-40425 red	
STRAIGHT - 90 CM		STRAIGHT - 76 CM	STRAIGHT - 66 CM
			
250 kg Item no.: 40-40402 grey Item no.: 40-40407 red	330 kg Item no.: 40-40430 grey Item no.: 40-40435 red	250 kg Item no.: 40-40401 grey Item no.: 40-40406 red	250 kg Item no.: 40-40400 grey Item no.: 40-40405 red


The anti-slip surface is available in grey (RAL 7037) and red (RAL 3003).

ACCESSORIES


Toilet roll holder
Item no.: 40-40452



Height adjustment
(Quick release)
Item no.: 40-40450



Height adjustment
(Adjusted with tool)
Item no.: 40-40451



NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS / ISO 17966:2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark – Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

ROPOX™
BETTER WAYS TO BETTER DAYS

– A PART OF AddLife®