



Ropox Mobilio Changing Bed

Ropox Mobilio Changing Bed is mobile and with brake wheels. This ensures that the changing bed can be fixed and is easy to move if there is a need for this.

Ropox Mobilio Changing Bed is connected to an electrical outlet with a three meter long cord for electric height adjustment. Height is infinitely adjustable from 65-100 cm using the hand control with spiral cord.

Ropox Mobilio Changing Bed are in two models with and without side rails as needed.

It is possible to buy mattresses for the table which provides added comfort and safety. The mattress has a clean surface, making it easy to clean after use.

There can be mounted baskets on rails under the table top for storage.

TECHNICAL INFORMATION

Ropox Mobilio Changing Bed

Lenght	Width 70	Width 80
90	40-30401	40-30419
100	40-30402	40-30420
110	40-30403	40-30421
120	40-30404	40-30422
130	40-30405	40-30423
140	40-30406	40-30424
150	40-30407	40-30425
160	40-30408	40-30426
170	40-30409	40-30427
180	40-30410	40-30428
190	40-30411	40-30429

All measuremets are in cm

Height: 65-100 cm

Max. weight load: 120 kg

Accessories:

Mattresses can be supplied in 3 sizes.

All mattresses are 3 cm and can be supplied with an non-slip surface.

Table top lenght x width	Mattress Item nr.
90 x 70	40-25793
100 x 70	40-25766
110 x 70	40-25768
120 x 70	40-25760
130 x 70	40-25761
140 x 70	40-25795
150 x 70	40-25762
160 x 70	40-25796
170 x 70	40-25763
180 x 70	40-25764
190 x 70	40-25765

All measuremets are in cm

Table top lenght x width	Mattress Item nr.
90 x 80	40-25769
100 x 80	40-25770
110 x 80	40-25771
120 x 80	40-25772
130 x 80	40-25773
140 x 80	40-25774
150 x 80	40-25775
160 x 80	40-25776
170 x 80	40-25777
180 x 80	40-25778
190 x 80	40-25779

Bed Guard -end wall
56 cm solid
V-nr. 40-25709

Bed Guard -front
70 cm for length 90-120 cm
V-nr. 40-25970

Bed Guard -front
80 cm for length 90-120 cm
V-nr. 40-25971

Bed Guard -front
70cm for length 120-190 cm
V-nr. 40-25972

Bed Guard -front
80 cm for length 120-190 cm
V-nr. 40-25973

Basket
1 pcs.
V-nr. 40-25714

Basket
2 pcs.
V-nr. 40-25715

Battery
Battery package incl.
wire intention
V-nr. 40-30762

Central lock
Surcharge for length 90 to 190
V-nr. 40-30750 to 40-30760

ABOUT US

Ropox develops, manufactures and markets products for installation in the homes of physically disabled people.

Through many years of collaboration with users and treatment groups we have acquired expertise in the design of solutions and products for:

- Shower rooms
- Kitchens
- Tables
- Nursing tables

To get more information visit our website www.ropox.com

No.1962164 1215

ROPOX A/S
RINGSTEDGADE 221
4700 NÆSTVED
TEL.: +45 5575 0500
INFO@ROPOX.DK
WWW.ROPOX.COM

ROPOXTM
BETTER WAYS TO BETTER DAYS