

”The optimum washbasin for both seated and standing users”



ROPOX SLIMLINE

WHEN YOU PREFER WASHBASIN AND MIRROR IN A COMBINED, SLIM AND ELEGANT SOLUTION

The SlimLine washbasin combines mirror and basin in one unit, and the mirror is therefore always at the right height no matter if you are sitting in a wheelchair or standing up. You adjust the height with a push on the control switch which is always within reach as it is located right under the mirror.

There are three types of washbasins to choose from, and each is particularly suited for specific users' needs:

- Standard: A simple and stylish washbasin for users without any special needs
- Support: A washbasin with good access and integrated handles that increases safety for users with limited physical capability
- Hospital: Porcelain washbasin with an inwardly curved front for hospitals

If a user's needs change over time, it is easy to replace the washbasin with another model since they are quite straightforward to install. The round corners of the SlimLine models minimize bacteria growth allowing you to keep a high hygiene level.



ROPOX SLIMLINE



Height adjustment:

70-100 cm

Max. user weight:

200 kg
(DS/EN17966-2016)

VARIANT/MODEL



Available with or without wash basin.

Mixer tap and plumbing kit as well as cover for installation is included with the washbasins.

ACCESSORIES




Plumbing kit
Item no.: 40-15612

Fixed light length 600mm
Item no.: 40-44182

Cover
Item no.: 40-15610

Hand control
Item no.: 40-15630

Plug with overflow
Item no.: 40-40935



NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: DS/ISO 17966-2016

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

ROPOXTM
BETTER WAYS TO BETTER DAYS

No. 1962238 0218