

”Doing more personal care
without assistance ”



ROPOX TOILET LIFTER

GET ON AND OFF THE TOILET MORE EASILY AND SAFELY

Most people would like to take care of their personal hygiene without assistance from others. The toilet lifter is one of the aids that can enable you to become more self-sufficient in getting on and off the toilet because it makes it easier and more effortless to sit down and stand up.

The toilet lifter lowers and lifts the user to the best sitting or standing position, and it is straightforward to set to the optimum height for both user and helper.

The round edges are gentle to the skin, and you do not get scratches, and the cleaning-friendly design makes the toilet lifter easy to keep clean ensuring a high level of hygiene.






ROPOX TOILET LIFTER

Height adjustment: Electric / Manual 40 cm
Basic 20 cm

Maximum load: 600 kg (max user load 400 kg
cf. DS/EN 17966:2016)

VARIANT/MODEL


TOILET LIFTER		
		
Electric Item no.: 40-45020	Manual Item no.: 40-45030	Basic Item no.: 40-45040

Toilet lifter unit in bottom position: Model Electric/Manual 118 x 79 x 20,5 cm and Model Basic 118 x 79 x 16,5 cm.

Dual flushing system at 3 and 6 litres.

ACCESSORIES


Toilet supports with height adjustable unit. See data sheet for arm supports for more options.




For optimum choice of toilet, please contact Ropox.



Upgrade kit (Manual to Electric)
Item no.: 40-45102



Long backrest
Item no.: 40-45100



Short backrest
Item no.: 40-45101



NB! Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: EN/EN17966-2016.

ROPOX A/S · Ringstedgade 221 · 4700 Næstved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

ROPOXTM
BETTER WAYS TO BETTER DAYS

No. 1962244 0218