

"More independence  
with your personal care"



## ROPOX TOILET SUPPORT ARMS

### GET SOMETHING TO LEAN ON

Toilet support arms are a great aid if you need to get on and off the toilet safely. Their ergonomic shapes provide a comfortable grip, and the cleaning-friendly surfaces optimize hygiene.

Because they are designed with two heights and an anti-slip surface, you have a firm grip from both standing and seated position. If you need a little extra support for getting up and sitting down, you can mount the standing/sitting attachments.

The support arms can be adjusted in height and at different angles to suit the individual user's needs.

If you need extra room in the bathroom, the support arms can be folded up to vertical position.







# ROPOX TOILET SUPPORT ARMS

Height adjustment: 12 cm




Maximum load: 400 kg (see max. user load below cf. DS/EN 66:2016)

## VARIANT/MODEL

WAVE - 90 CM		
		
80 kg Item no.: 40-40115	135 kg Item no.: 40-40125	190 kg Item no.: 40-40145

WAVE - 76 CM

110 kg Item no.: 40-40151

Can also be used for Ropox Toilet Lifter

STRAIGHT - 90 CM		
		
80 kg Item no.: 40-40110	135 kg Item no.: 40-40120	190 kg Item no.: 40-40140


STRAIGHT - 76 CM

110 kg/190 kg Item no.: 40-40150/40-40152


Can also be used for Ropox Toilet Lifter

## ACCESSORIES

Height adjustment  
(adjusted with tool)  
Item no.: 40-40906



Height adjustment  
(adjusted without tool)  
Item no.: 40-40907



Toilet roll holder  
Item no.: 40-40910



Standing/Sitting attachment  
Item no.: 40-40925



**NB!** Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.

The items are produced according to the following Standard: EN/EN17966-2016.

ROPOX A/S · Ringstedgade 221 · 4700 Naestved, Denmark · Phone: +45 5575 0500 · info@ropox.com · www.ropox.com

**ROPOX**<sup>TM</sup>  
BETTER WAYS TO BETTER DAYS

No. 1962243 0218